


# Test Report

- **Certificate No. :** TE-20-04011
- **Client**  
Corporation Name :CSM IMPLANT  
Address : B205, Techno B/D, Kyung-pook National Univ. 47, Gyeongdae-ro 17gil
- **Sample Description:** C5F4414 AH0410-1
- **Date of Receipt :** 2020. 12.08
- **Date of Test :** 2020. 12.15
- **Test method used :** Reference the "Test method" in the next page(s)
- **Test Results :** Reference the "Test Results" in the next page(s)

**Tested by**  
**Name :** PARK JIN  
 HONG



**Technical Manager**  
**Name :** Son young bum



· The above test report is the result of the sample which a client submitted and should be used only it's own object

2021.08.18



**President, Daegu Mechatronics & Materials Institute**




# Test Results

Certificate No. : TE-20-04011

( 2 ) Page of ( 2 ) Pages



## 1) Machine information

Test Machine	Model	Test Area	Test Method
Noncontact Roughness Test Machine	NT-8000	(0.208 x 0.156) mm	Noncontact Roughness Test

## 2) Test Result

Item Name	Test Item	Unit	Test Result
C5F4414 AH0410-1	Roughness (Ra)	$\mu\text{m}$	1.88
			1.72
			1.90

End.

