


Test Report

- **Certificate No. :** TE-20-04012
- **Page No. :** 1 Page of 2
- **Client**
Corporation Name :CSM IMPLANT
Address : B205, Techno B/D, Kyung-pook National Univ. 47, Gyeongdae-ro 17gil
- **Sample Description:** LC5F4810 AH1708-1
- **Date of Receipt :** 2020. 12.08
- **Date of Test :** 2020. 12.15
- **Test method used :** Reference the "Test method" in the next page(s)
- **Test Results :** Reference the "Test Results" in the next page(s)

Tested by
Name : PARK JIN
 HONG



Technical Manager
Name : Son young bum



· The above test report is the result of the sample which a client submitted and should be used only it's own object

2021.08.18



President, Daegu Mechatronics & Materials Institute




Test Results

Certificate No. : TE-20-04012

(2) Page of (2) Pages



1) Machine information

Test Machine	Model	Test Area	Test Method
Noncontact Roughness Test Machine	NT-8000	(0.208 x 0.156) mm	Noncontact Roughness Test

2) Test Result

Item Name	Test Item	Unit	Test Result
LC5F4810 AH1708-1	Roughness (Ra)	μm	10.64
			11.83
			11.83

End.

